Recent Advances in Application of Cerebral Oximetry in Adult Cardiovascular Surgery

Seminars in Cardiothoracic and Vascular Anesthesia Volume XX Number X Month XXXX xx-xx © 2008 Sage Publications 10.1177/1089253208316443 http://scv.sagepub.com hosted at http://online.sagepub.com

Gregory W. Fischer, MD

Cerebral oximetry is a noninvasive technology that continuously monitors cerebral tissue oxygen saturation, which is a sensitive index of global cerebral hypoperfusion. On the basis of near-infrared spectroscopy technology, information is provided on the availability of oxygen in brain tissue at risk during numerous pathological conditions. Complementary to the arterial oxygen saturation measured by pulse oximetry, cerebral tissue oxygen saturation reflects regional cerebral metabolism and the balance of local cerebral oxygen supply/demand. Recently, patient management with the guidance of cerebral oximetry has resulted in improved patient outcomes. This review will briefly describe the physics behind cerebral oximetry and will provide an overview of the literature focusing primarily on articles published within the past 5 years.

Keywords: NIRS; cardiac surgery; cerebral oximetry

erebral oximetry, based on near-infrared spectroscopy (NIRS) technology, provides information on the availability of oxygen in brain tissue at risk during numerous pathological conditions.¹ Cerebral oximetry measures regional cerebral tissue oxygen saturation (SctO₂) at the microvascular level (arterioles, venules, and capillaries only).²⁻⁴ Complementary to the arterial oxygen saturation (SaO₂) measured by pulse oximetry, SctO₂ reflects regional cerebral metabolism and the balance of local cerebral oxygen supply/demand. The advantages of cerebral oximetry are: (1) it provides SctO₂ values continuously and noninvasively at the bedside,⁵ and (2) SctO₂ is a sensitive index of cerebral hypoxia and/or cerebral ischemia,6 which is one of the main causes of brain injury in clinical settings.7

Cerebral oximetry, otherwise know as NIRS, is the subject of many reviews.^{1,5,8} In one of the most recent papers, Taillefer and Denault⁹ systematically reviewed the clinical efficacy of cerebral oximetry monitoring in adult cardiac surgery. On the basis of reviews of publications dated pre-January 2004, they concluded that further investigation is needed to prove the clinical utility of NIRS in cardiac surgery. In the past few years, there have been new advances in NIRS technology and in clinical studies of cerebral oximetry in adult cardiovascular surgery. This review will focus on these recent developments and our experiences with cerebral oximetry monitoring during adult cardiovascular surgery at our institution.

The Significance of Cerebral Oximetry Monitoring for Cardiovascular Surgery

Monitoring brain oxygenation is critical in providing information to guide patient management in many clinical situations.^{10,11} Currently, brain oxygenation can be measured invasively by jugular bulb oximetry $(SjvO_2)$ or brain tissue Po₂ sensor.^{12,13} However, these modalities are invasive and difficult to use. Therefore, a bedside cerebral oximeter that can provide noninvasive measurement of cerebral oxygenation is highly desirable.

Despite a decline in the overall mortality after coronary artery bypass grafting (CABG) and valvular surgery with cardiopulmonary bypass (CPB), the rates of cognitive dysfunction have not improved.¹⁴ In some reports, most patients had subtle signs of impaired

From the Department of Anesthesiology, Mount Sinai Medical Center, New York, New York.

Address correspondence to: Gregory W. Fischer, MD, Department of Anesthesiology, Mount Sinai Medical Center, PO Box 1010, One Gustave L Levy Place, New York, NY 10029; e-mail: Gregory .Fischer@mountsinai.org.

cognitive performance, with incidences ranging from 60% to 80%.15 There are 2 different forms of brain injury that may occur after CABG or thoracic aorta surgery: neurological dysfunction (ND) and neurocognitive dysfunction (NCD). The ND is defined as clinically evident focal or global neurological injury resulting in stroke, hypoxic encephalopathy, transient ischemic attack, or stupor. The NCD is defined as postoperative confusion, agitation, delirium, prolonged obtundation, or transient Parkinsonism. The NCD occurs more frequently, affecting 40% to 80% of CABG or aortic surgery patients, depending on the method of detection used.¹⁶ Although it is easy to diagnose postoperative ND, NCD is more subtle and needs to be evaluated with a full battery of neurocognitive testing administered by trained professionals. Recent reports based on accurate neurocognitive testing before and after surgery have suggested that NCD can no longer be considered a benign self-limiting condition but rather a long-lasting neurocognitive insult capable of reducing quality of life by impairing memory and fine motor function.¹⁷

The etiology of brain injury following cardiac surgery is still not completely understood and is somewhat controversial. Some of the possible mechanisms include diffuse microembolization, cerebral hypoperfusion, and metabolic factors. The incidence of injury seems to be higher when CPB duration exceeds 70 minutes and when there is rapid rewarming, particularly in the older surgical population.¹⁸ Regardless of the immediate cause, such persistent cognitive dysfunction likely results from brain ischemia during surgery, which may be a result of focal arterial embolism,¹⁹ global hypoperfusion of the brain,²⁰ or an interaction of the 2.²¹

Strategies for preventing arterial embolism and brain hypoperfusion differ. To avoid arterial embolism during CPB, arterial line filters, intraoperative imaging, and careful manipulation by the surgeon are essential. However, diffuse hypoperfusion of the brain can be avoided only by very careful planning; if something goes wrong, the only hope is early detection and immediate restoration of adequate perfusion before irreversible brain damage develops. For this purpose, sensitive, real-time monitoring of brain ischemia during such surgical procedures is needed.²² At the present time, cerebral oximetry can fulfill this role as the only feasible technology that monitors cerebral hypoxia and/or cerebral ischemia noninvasively and continuously.

The Evolution of Cerebral Oximetry

The idea of using NIRS to measure cerebral tissue oxygenation noninvasively was first introduced by Jobsis in 1977.²³ The principle of NIRS is based on the known fact that near-infrared light passes through skin and skull readily and is absorbed by certain biological molecules in the brain.^{24,25} The NIRS measures oxygenation in biological tissue mainly at the microcirculation level (capillaries, arterioles, and venules) based on different absorption characteristics of the chromophores oxyhemoglobin (Hbo₂) and deoxyhemoglobin (Hb) in the near-infrared spectrum. A "biological spectroscopic window" exists at the wavelength range 660 to 940 nm because only a few chromophores like Hb and Hbo, strongly absorb light in this spectra range, allowing light to penetrate tissue to a great distance. In this wavelength range, absorption of light due to other biological compounds and tissues, such as water, lipids, skin, and bone, is lower in magnitude, and these biological compounds generally have a flat absorption spectra when compared with Hb and Hbo₂.

Over the years, many NIRS devices (cerebral oximeters) were developed by either research groups or commercial entities.²⁶⁻³⁰ Early versions of NIRS monitors have limited capacity or suffer from poor performance. With the advances in optical and electronic technologies and by the improved understanding of the underlying basic theory and model of light propagation in biological tissue, NIRS has evolved to semi-quantitative devices⁸ and finally to absolute quantitative instruments.

Technically, NIRS instrumentation implements continuous wave (CW), time domain (TD), or frequency domain (FD) technology^{31,32} often combined with spatially resolved (SR)33 methods. Most CW NIRS instruments measure changes in light absorption by recording the changes in the intensity of light both transilluminated through or reflected by photon scattering in tissue. The CW instruments usually do not measure the path length. Instead, path length is canceled out mathematically in oxygen saturation algorithms or is substituted with a predetermined average path length obtained from data collected by TD or FD instruments. The TDtype instruments measure the temporal response of tissue to an ultrashort input light pulse that lasts only for a few picoseconds. The temporal distribution of light emerging from the tissue surface is then

detected with a streak camera or photon counter.^{34,35} The flight times relative to a reference pulse of individual photons are measured to obtain a histogram of the distribution of arrival times (temporal point spread function [TPSF]). By fitting the TPSF to a light transport model, values of absorption and scattering (μ_a and μ_s) can be calculated.³⁶ The TD-type NIRS instruments are considered to be more accurate, but they are large, expensive, and confined to dedicated optical laboratories. In the FD-type instrument, the light source is intensity modulated at radio frequencies (100~200 MHz). For light exiting the tissue, the light intensity I, phase shift Φ , and modulation depth M relative to the input signal are measured. By fitting the data to a diffusion model, the changes in chromophore concentrations can be calculated from attenuation changes.³⁷ The SR refers to a methodology used in some NIRS devices by measuring intensity changes at multiple light sources and detector spacings. The SR devices use the gradients of the optical signal from this multidistance measurement to quantify the NIRS parameters.³⁷

Regardless of the type of methods applied, a NIRS instrument comprises the means to deliver infrared light to the object under study and the detecting means to collect the light after it propagates through the object. For a clinical NIRS monitor, the challenges center around the key components, which are (a) the optical probe that is easily attached to the object and is able to safely deliver energy to the object under study, as well as to collect the very low-intensity energy after it propagates through the object; and most of all, (b) the algorithm that processes the received signals and calculates the NIRS parameters.

At present, there are only a few cerebral oximeters commercially available. Of these, only 2 companies have food and drug administration (FDA) clearance to market cerebral oximeters in the United States.

INVOS (Somanetics Corporation, Troy, Michigan)

The INVOS cerebral oximeter has been commercially available for over 10 years.

The INVOS uses 2 wavelengths of infrared light (730 and 810 nm) from light-emitting diodes (LEDs).³⁸ This device has 2 channels for bilateral brain monitoring. The INVOS has disposable sensors, which contain LEDs and 2 light detectors at

fixed distances from the light source (3 and 4 cm for adult sensor). The INVOS uses an SR algorithm assuming that there is no wavelength dependence of scattering and that there is linearity over the 1 cm between the 2 detectors.³⁹ The device calculates and displays the value of regional cerebral oxygen saturation (rSO₂). According to the FDA "Indications for Use," the INVOS only provides trend monitoring of regional hemoglobin oxygen saturation of blood in the brain of an individual.²⁹ This means that clinical intervention should be based on changes of rSO₂ from the initial baseline value, which is usually obtained at preinduction of anesthesia. In other words, the value of rSO₂ from the INVOS cerebral oximeter should be viewed with caution if the preinduction value is unknown. In fact, because we use cerebral oximetry routinely in cardiac surgery cases, the author has personally observed events during which the rSO₂ value of the INVOS was in the thirties (which is consider gravely low⁴⁰ if taken as an absolute value) in alert patients at preinduction.

Trend-only cerebral oximetry may be inadequate. Some studies suggest that perioperative management can be modified based on the application of trend-only cerebral oximetry. Because these monitors measure trends only, a baseline has to be established first, and cerebral oximetry values need to be maintained at or near the preoperative baseline.^{41,42} Another approach is to keep the cerebral oxygen saturation at levels within 20% to 25% of the anesthesia preinduction value.43,44 However, studies have shown that 75% of patients undergoing CABG have a significant impairment in baseline regional cerebral perfusion (rCP).45 Other studies also demonstrated that patients undergoing cardiovascular surgery have a high prevalence of cerebral vascular disease in varying degrees.⁴⁶ In addition, abnormal preoperative rCP was found to be a strong indicator for postsurgical decline in neuropsychological testing.⁴⁷ These findings suggest that it is difficult to define a normal preinduction baseline value for the trendonly cerebral oximetry. A percentage drop based on the unreliable baseline value is more questionable to serve as a threshold for clinical intervention. In fact, studies have confirmed that although a trend-only cerebral oximeter can detect adverse brain oxygenation by measuring the change in rSO₂ from a baseline value, it cannot provide accurate, reliable, normal, and threshold values of SctO₂.⁴⁸⁻⁵⁴

FORE-SIGHT (CAS Medical Systems Inc, Branford, Connecticut)

The FORE-SIGHT absolute cerebral oximeter is a new device that has recently emerged in the market.⁵⁵The FORE-SIGHT monitor was developed with research grants from the National Institute of Neurological Disorders and Stroke of the National Institute of Health (NIH).⁵⁶ The FORE-SIGHT monitor is a continuous wave, spatially resolved, near-infrared spectrometer that measures absolute cerebral tissue oxygen saturation (SctO₂%). This device has 2 channels for bilateral brain monitoring. The disposable sensor has a light source composed of a fiber optic light guided from the monitor with a prism and 2 light detectors (scalp and brain) at fixed distances (1.5 and 5 cm) from the light source. The scalp detector samples returning light from the extracranial tissue, whereas the brain detector samples returning light signal from both the brain and the extracranial tissue. The signal from the scalp detector is used to cancel extracerebral interference from the signal of the brain detector to obtain information that mostly comes from the brain. Laser light at an FDAdeemed safe level (Class I laser) is projected into the brain in 4 precise (bandwidth < 1 nm) wavelengths (690, 780, 805, and 850 nm) to capture information needed for the algorithm to calculate the absolute value of SctO₂. From the NIH funded research,⁵⁶ it was determined that 4 precise wavelengths are needed to maximize the measurement accuracy of oxyhemoglobin and deoxyhemoglobin in determining absolute SctO₂, by compensating for wavelength-dependent scattering losses and by accounting for interference from other background light absorbers⁵⁷ (such as fluid, tissue, and skin pigmentation). Because the FORE-SIGHT provides an absolute value of SctO₂, it is possible to establish threshold values for SctO₂ to guide clinical interventions. Furthermore, no preinduction baseline measurement is needed for obtaining meaningful readings.

FORE-SIGHT absolute cerebral oximeter determined SctO₂ is defined as the ratio of concentrations of Hbo₂ to Hb + Hbo₂ in the brain tissue. The value of SctO₂ reflects a proportional mix of arterial and venous blood, which was validated from arterial and internal jugular venous blood.⁵⁸ It is estimated that the NIRS cerebral oximeter interrogated brain tissue microvasculature is approximately 70% venous and 30% arterial during most physiological conditions in humans based on the positron emission tomographic studies.⁵⁹ In validation studies, FORE-SIGHT–determined $SctO_2$ showed a strong correlation with invasive blood samples over the spectrum of pulse oximeter–determined arterial oxygen saturation (SPo₂) values between 70% and 100% from healthy volunteers.^{60,61} The bias and precision (standard deviation [SD] = 1) for the FORE-SIGHT SctO₂ compared with reference SctO₂ derived from co-oximetry of arterial and jugular bulb blood samples was 0.18 ± 3.7 (SD = 1). It is clinically accepted that SjvO₂ has a normal lower safe limit of ~45%.^{62,63} The FORE-SIGHT SctO₂ over a wide range of oxygen saturation values. Therefore, the absolute FORE-SIGHT cerebral oximeter lower safe SctO₂ threshold is estimated to be approximately 55%.

MacLeod et al⁶⁴ published an abstract containing the initial results of the FORE-SIGHT study on cardiac patients. In all, 33 patients presenting for CABG and/or valvular surgery were included. The median SctO₂ value in patients prior to induction was 70% with small variability (SD = 4.1%). This preinduciton SctO₂ value was independent of age, race, or sex. After initiation of CPB, a decrease in SctO₂ was seen (65 ± 5%), but eventually SctO₂ recovered to close to preinduction values by the time of chest closure (68 ± 5%). Interestingly, patients undergoing valve procedures had lower SctO₂ values for longer duration while on CPB.

Fischer et al⁶⁵ presented preliminary results of the use of the FORE-SIGHT as a monitor of cerebral oxygenation during aortic surgery with deep hypothermic circulatory arrest. The postinduction SctO₂ was $68.5 \pm$ 6%. After the initiation of CPB and cooling of the patient, the SctO₂ value increased in most patients (84% ± 4.9%). The rate of SctO₂ increase matched the rate of core cooling. After initiation of circulatory arrest, a decline in SctO₂ was seen (59 ± 8.9%). After reestablishing cerebral perfusion, the SctO₂ values quickly returned to prearrest levels.

Applications of Cerebral Oximetry Monitoring

Two of the main areas of interest for the clinician are the ability of the device to predict a catastrophic event and to have a direct positive effect on clinical outcome. Both of these goals are ambitious, and the outcome goal is particularly difficult to prove. However, studies are emerging revealing that optimization of SctO₂ improves overall morbidity and decreases hospital length of stay. The latter endpoint is appealing not only to patients and their physicians but also to hospital administrators and society in general.

Catastrophic Events Avoidance

The discipline of anesthesiology is frequently compared with that of aviation. Both fields are historically viewed by the layman as extremely dangerous. However, with advances in science and safety algorithms occurring over the course of the past decades, both anesthesiology and aviation have become extremely safe. Unfortunately, despite rigorous training and numerous checklists, unexpected catastrophic events still occur. Cerebral oximetry adds an additional safety check for the clinician during the perioperative period. Catastrophic events that may occur during the course of surgeries, which would otherwise have gone unrecognized resulting in permanent damage or death, can now possibly be identified early enough for corrective measures to be implemented. The author will review the literature and describe his own experiences regarding the catastrophic events that could be avoided by the continuous monitoring of $SctO_2$.

Most reports regarding catastrophic events are case studies focused around the pediatric population. Rodriguez et al66 and Sung and Chong67 reported on the utility of cerebral oximetry to detect superior vena cava obstruction caused by malpositioned venous cannulae. In addition to cerebral oximetry, Rodriguez et al⁶⁶ used transcranial Doppler and electroencephalographic (EEG) monitoring. During upper venous congestion, all monitoring modalities revealed pathological values. However, the report of Sung and Chong⁶⁷ relied solely on cerebral oximetry to detect venous cannula malpositioning. A severe decline reading cerebral oximetry and plethora of the face were the only 2 signs of an ensuing catastrophic event. All other physiological values (eg, blood pressure, Pao₂, and Paco₂) were normal.

Yeh and Austin⁶⁸ reported on the use of a multimonitoring approach to detect catastrophic neurological events. During the fenestration phase of a Fontan procedure, air was accidentally entrained into the systemic circulation. An acute deterioration in the readings of transcranial Doppler, cerebral oximetry, and EEG monitoring was detected. Corrective interventions were swiftly implemented. No neurological sequelae resulted from the air embolism. Orihashi et al⁶⁹ report on the utility of cerebral oximetry to detect malpositioned cerebral perfusion catheters used during selective antegrade cerebral perfusion for aortic arch repair. In 4 of 35 consecutive cases, a malpositioned catheter could be detected by an unexplained decline in cerebral oximetry.

Similar use of cerebral oximetry regarding correct placement of cerebral perfusion catheters is mentioned by Hagl et al⁷⁰ in a detailed review of different perfusion and monitoring techniques used during hypothermic circulatory arrest.

At the author's institution, 2 catastrophic events were first detected by the FORE-SIGHT absolute cerebral oximeter. A patient presenting for aortic arch surgery revealed an unusual SctO₂ pattern shortly after CPB was initiated. A sharp decline in $SctO_2$ (40%) was observed despite cooling to 20°C. Echocardiography monitoring indicated that the low SctO₂ value seen in this patient might have been caused by air emboli. The surgeon was informed, and the field was submerged with saline. The source of air entrainment was believed to originate from a left ventricular vent entering the upper right pulmonary vein. The insertion site was snared down with a purse-string suture. Subsequent air entrainment could not be identified thereafter. This patient was extubated on postoperative day 5 (compared with an average extubation time of 1 or 2 days) due to severe confusion and lack of orientation to time, place, or person. Long-term neurological follow-up evaluation for all patients is pending.

A second case in which cerebral oximetry gave the clinician an early warning sign was during an elective Bentall procedure. During the rewarming phase, after completion of the procedure, a sudden unexpected drop in SctO₂ was detected (SctO₂ 53%, Figure 1). An arterial blood gas analysis revealed respiratory alkalosis (Paco₂ 26 mm Hg). The cerebral desaturation was due to hypocapnia, which led to cerebral vasoconstriction and hypoperfusion. The event occurred because mechanical ventilation was commenced without decreasing the fresh gas delivery to the membrane oxygenator of the cardiopulmonary bypass circuit while still on partial CPB. As soon as it was recognized and the respiratory derailment corrected, SctO₂ promptly returned to the normal values. No neurological sequelae resulted from this intraoperative event.



Figure 1. An event of cerebral desaturation due to hypocapnia was quickly recognized with cerebral oximetry monitoring (FORE-SIGHT).

The event occurred because mechanical ventilation was initiated without decreasing the fresh gas delivery to the membrane oxygenator of the cardiopulmonary bypass circuit while still on partial cardiopulmonary bypass. SctO₂ indicates cerebral tissue oxygen saturation; DHCA, deep hypothermic circulatory arrest; CPB, cardiopulmonary bypass.

Tailored Patient Management

An increasing number of evidence is mandating the need for a tailored patient management protocol. Current approaches for managing pump flow, arterial blood pressure, and pH during the cardiac surgery are based on studies that included a few elderly or high-risk patients and predated many other contemporary practices.⁷¹ For example, watersheddistribution stroke happens more frequently in patients undergoing cardiac surgery than in the general stroke population (over 40% vs 2%-5%, respectively).⁷²⁻⁷⁴ Gottesman et al⁷⁵ reported that the mechanism of watershed stroke after cardiac surgery may include an intraoperative decline in blood pressure from a patient's baseline. This suggests that following a standard protocol to maintain blood pressure at a predetermined value during cardiac surgery is insufficient for some patients. Real-time cerebral oximetry is capable of evaluating the effects of blood pressure and perfusion changes to cerebral oxygenation. Therefore, cerebral oximetry can be used for this tailored patient management approach.

A growing body of literature is emerging that demonstrates improved outcomes in patients in whom cerebral tissue oxygenation was monitored and optimized. In a large, nonrandomized study of patients undergoing cardiac surgery, Goldman⁴¹ reported that the treatment group using cerebral oximetry monitoring to maintain regional cerebral oxygen saturation (rSO₂, INVOS, Somanetics) values at or near the patient's preinduction baseline demonstrated a significantly lower incidence of permanent stroke⁶⁷. The weakness of this study is its nonrandomized and retrospective design. The control group consisted of a historical cohort. Patients who underwent cardiac surgery 18 months prior to the introduction of cerebral oximetry at the author's institution were subsequently compared with the treatment group consisting of patients who underwent cardiac surgery the following 18 months. Additionally, the author acknowledges that they were unable to determine how close rSO₂ values were maintained to baseline values during the course of surgery. Despite these shortcomings, one could argue, as the author pointed out, that better compliance with a predefined target value would translate into potentially better outcomes. The fact that the study group consisted of patients with more comorbidities, yet had a lower incidence of permanent stroke; decreased need for prolonged ventilation; and shorter length of hospital stay is impressive.

Murkin et al⁷⁶ published an article in January 2007 that corrected many of the shortcomings of the Goldman⁴¹ paper. A randomized, blinded study of 200 CABG patients demonstrated prospectively that rSO₂ (INVOS, Somanetics) monitoring is associated with a significant improvement in overall outcome after cardiac surgery. Murkin utilized an intraoperative management protocol designed to maintain rSO₂ values at or above 75% of the preinduction baseline. The results were associated with a significant improvement in overall organ function and decreased postoperative length of stay. Fewer patients suffered a stroke in the intervention group compared with patients in the control group. However, a significant reduction in stroke rate could not be seen due to underpowering of the study. Nevertheless, the fundamental hypothesis that cerebral oxygenation could be used as a surrogate marker for systemic perfusion and improvement in clinical outcome was supported by the results of this study.

Orihashi et al⁷⁷ showed in an observational study examining 59 consecutive patients presenting for aortic surgery using selective cerebral perfusion (SCP) that a sustained decrease in rSO_2 (TOS-96, Tostec CO, Tokyo, Japan) was observed in patients with infarcts, suggesting hypoperfusion as the etiology. Transient neurological events occurred more frequently in patients with sustained drops below 55% for over 5 minutes. Recommendations were made that a decline in rSO_2 below 55% should be addressed without delay.

In another randomized, prospective study of elderly patients undergoing major abdominal surgery, Casati et al⁷⁸ reported on improved mini mental scores (MMS), more rapid post anesthesia care unit discharge, and shortened length of hospital stay in 122 geriatric patients after major abdominal surgery in which cerebral oximetry was monitored and optimized. Although the utilization of the MMS to address postoperative cognitive impairment is crude, this study did demonstrate that 20% of all patients experienced a decrease in rSO₂ (INVOS, Somanetics) below 75% of baseline without any change in systemic oxygenation (SPo₂).

Our experience with the FORE-SIGHT absolute cerebral oximeter also demonstrated that this device can be used as a guide for patient management during cardiac surgery. Fischer et al⁷⁹ reported an unexpected SctO₂ pattern in a patient presenting for repair of an ascending aortic aneurysm. Despite uneventful hemodynamic management (mean arterial pressure [MAP] 50 mm Hg) on CPB a steady decline in SctO₂ was observed. Postoperative analysis showed that unlike other patients who have normal patterns of SctO₂ (Figure 2), in this patient, $SctO_2$ variation corresponded to dynamic fluctuations of mean arterial pressure MAP (Figure 3). This suggests that the failure of SctO₂ to rise with cooling was attributed to the loss of cerebral autoregulation. These findings indicate that hemodynamic management on CPB needs to be tailored to the individual patient. Maintaining MAP at 50 to 60 mm Hg during hypothermic CPB is tolerated by most patients, but this level seems to be inadequate for certain patients.

However, not all studies have shown a positive correlation between cerebral oximetry monitoring and postoperative outcomes. Reents et al⁸⁰ report in an observational study of 47 patients undergoing CABG that no correlation could be found between intraoperative regional rSO₂ values as assessed by the INVOS 4100 device and early postoperative cognitive performance. Careful review of this study does alert the reviewer to flaws in its design. The author uses 2 cutoff points as markers for low cerebral perfusion, a decrease in rSO₂ values under 40% and/or a decrease of more than 25% from baseline values. The INVOS is approved as a trend-only cerebral oximeter. Therefore, expecting this device to identify an absolute value, which would represent a



Figure 2. Normal pattern of $SctO_2$ of a patient undergoing aortic surgery with DHCA.

The SctO₂ increased with cooling to 20°C and increased further with hard cooling to 13°C. Fluctuation of MAP, which was maintained around 50 mm Hg, had no effect on SctO₂. SctO₂ indicates cerebral tissue oxygen saturation; DHCA, deep hypothermic circulatory arrest; MAP, mean arterial pressure; CPB, cardiopulmonary bypass.



Figure 3. Observation of a steady decline in $SctO_2$ with cooling despite uneventful hemodynamic management (MAP 50 mm Hg) on cardiopulmonary bypass.

The SctO₂ variation corresponded to dynamic fluctuations of mean arterial pressure MAP, suggesting that the failure of SctO₂ to increase with cooling was attributed to the loss of cerebral autoregulation. SctO₂ indicates cerebral tissue oxygen saturation; MAP, mean arterial pressure; CPB, cardiopulmonary bypass; DHCA, deep hypothermic circulatory arrest.

threshold correlating to poor postoperative outcomes, must be viewed as a methodological flaw. Additionally, the author defined baseline value as the value obtained after induction of general anesthesia and transfer of the patient from the induction room to the operating theater. A true baseline value should be obtained prior to the induction of anesthesia as was done in the Murkin⁷⁶ paper.

Summary

Cerebral oximetry has emerged as an exciting new technology that many cardiac anesthesiologists have incorporated into their monitoring standard of care. In the past, critics were quick to point out the lack of evidence positively linking the use of cerebral oximetry with patient outcomes. However, recent randomized prospective studies have confirmed the clinical efficacy of cerebral oximetry monitoring. As with any new technology, improvement in accuracy and reliability is an ongoing process. The past decade has seen tremendous advances in NIRS technology. The introduction of the FORE-SIGHT absolute cerebral oximeter in 2006 must be considered a groundbreaking achievement. Being able to measure absolute values of SctO₂ opens new possibilities for the clinician. Most importantly, it will now be possible to establish threshold values for SctO₂ to guide clinical interventions. Additionally, patients for whom a baseline value could not be established (intubated patients in emergency department or intensive care unit) can now profit from real-time cerebral oxygenation monitoring.

The past 5 years have seen the publication of well-designed studies clearly revealing that cerebral oxygenation, if monitored and optimized intraoperatively, can improve outcome and can shorten the hospital length of stay^{76,79}.

Although it is by no means the author's intent to endorse this technology solely on the basis of being an attractive concept, it must be stated that the overall risk to possible benefit clearly lies on the side of cerebral oximetry monitoring.

References

- Ferrari M, Mottola L, Quaresima V. Principles, techniques, and limitations of near infrared spectroscopy. *Can J Appl Physiol*. 2004;29:463-487.
- Benni PB, Chen B, Dykes FD, et al. Validation of the CAS neonatal NIRS system by monitoring VV-ECMO patients: preliminary results. *Adv Exp Med Biol.* 2005; 566:195-201.
- Kurth CD, Steven JL, Montenegro LM, et al. Cerebral oxygen saturation before congenital heart surgery. *Ann Thorac Surg.* 2001;72:187-192.

- Watzman HM, Kurth CD, Montenegro LM, et al. Arterial and venous contributions to near infrared cerebral oximetry. *Anesthesiology*. 2000;93:947-953.
- Madsen PL, Secher NH. Near-infrared oximetry of the brain. Prog Neurobiol. 1999;58:541-560.
- Yao FS, Tseng CC, Ho CY, Levin SK, Illner P. Cerebral oxygen desaturation is associated with early postoperative neuropsychological dysfunction in patients undergoing cardiac surgery. J Cardiothorac Vasc Anesth. 2004;18:552-558.
- Hogue CW Jr, Palin CA, Arrowsmith JE. Cardiopulmonary bypass management and neurologic outcomes: an evidencebased appraisal of current practices. *Anesth Analg.* 2006; 103:21-37.
- Owen-Reece H, Smith M, Elwell CE, Goldstone JC. Near infrared spectroscopy. *Br J Anaesth.* 1999;82:418-426.
- 9. Taillefer MC, Denault AY. Cerebral near-infrared spectroscopy in adult heart surgery: systematic review of its clinical efficacy. *Can J Anaesth.* 2005;52:79-87.
- Feldman Z, Robertson CS. Monitoring of cerebral hemodynamics with jugular bulb catheters. *Crit Car Clin.* 1997; 13:51-77.
- De Deyne C, Van Aken J, Decruyenaere J, Struys M, Colardyn F. Jugular bulb oximetry: review on a cerebral monitoring technique. *Acta Anaesthesiol Belg.* 1998;49: 21-31.
- al-Rawi PG, Hutchinson PJ, Gupta AK, Piechnik SK, Pickard JD, Kirkpatrick PJ. Multiparameter brain tissue monitoring—correlation between parameters and identification of CPP thresholds. *Zentralbl Neurochir*. 2000; 61:74-79.
- Jödicke A, Hübner F, Böker DK. Monitoring of brain tissue oxygenation during aneurysm surgery: prediction of procedure-related ischemic events. *J Neurosurg.* 2003; 98:515-523.
- Roach GW, Kanchuger M, Mangano CM, et al. Adverse cerebral outcomes after coronary bypass surgery. N Engl J Med. 1996;335:1857-1863.
- Newman MF, Kirchner JL, Phillips-Bute B, et al. Longitudinal assessment of neurocognitive function after coronary artery bypass surgery. N Engl J Med. 2001;344:395-402.
- Richard FM, Bill IW: Normothermic versus hypothermic cardiopulmonary bypass: central nervous system outcome. J Cardiothorac Vasc Anesth. 1996;10:45-53.
- 17. Newman MF, Grocott HP, Mathew JP, et al. Report of the substudy assessing the impact of neurocognitive function on quality of life 5 years after cardiac surgery. *Stroke*. 2001;32:2874-2881.
- Newman MF, Kirchner JL, Phillips-Bute B, et al. Longitudinal assessment of neurocognitive function after coronary artery bypass surgery. N Engl J Med. 2001;344: 395-402.
- Ergin MA, Galla JD, Lansman L, Quintana C, Bodian C, Griepp RB. Hypothermic circulatory arrest in operations on the thoracic aorta. Determinants of operative mortality

and neurologic outcome. J Thorac Cardiovasc Surg. 1994; 107:788.

- Plestis KA, Gold JP. Importance of blood pressure regulation in maintaining adequate tissue perfusion during cardiopulmonary bypass. *Semin Thorac Cardiovasc Surg.* 2001;13:170-175.
- 21. Hagl C, Khaladj N, Karck M, et al. Hypothermic circulatory arrest during ascending and aortic arch surgery: the theoretical impact of different cerebral perfusion techniques and other methods of cerebral protection. *Eur J Cardiothorac Surg.* 2003;24:371-378.
- 22. Murkin JM. Hemodynamic changes during cardiac manipulation in off-CPB surgery: relevance in brain perfusion. *Heart Surg Forum*. 2002;5:221-224.
- 23. Jobsis FF. Noninvasive, infrared monitoring of cerebral and myocardial oxygen sufficiency and circulatory parameters. *Science*. 1977;198:1264-1267.
- 24. Elwell CE. A Practical Users Guide to Near Infrared Spectroscopy. Hamamatsu, Japan: Photonics KK; 1995.
- McCormick PW, Stewart M, Ray P, Lewis G, Dujovny M, Ausman J. Neurological Res. 1991;13:65-70.
- 26. Cope M, Delpy DT. System for long-term measurement of cerebral blood and tissue oxygenation on newborn infants by near infra-red transillumination, *Med Biol Eng Comput.* 1988;26:289.
- McKeating EG, Monjardino JR, Signorini DF, Souter MJ, Andrews PJ. A comparison of the Invos 3100 and the Critikon 2020 near-infrared spectrophotometers as monitors of cerebral oxygenation. *Anaesthesia*. 1997;52: 136-140.
- Kurth CD, Thayer WS. A multiwavelength frequencydomain near-infrared cerebral oximeter, *Phys Med Biol*. 1999;44:727.
- 29. 510 (K) K001842, K960614, & K971628. Somanetics Invos 3100-5100 Adult/Pediatric Cerebral Oximeter and 5100 SSP and 5100 SPF Pediatric SomaSensor. Rockville, MD: FDA; 2005.
- Komiyama T, Onozuka A, Miyata T, Shigematsu H. Oxygen saturation measurement of calf muscle during exercise in intermittent claudication. *Eur J Vasc Endovasc Surg.* 2002;23:388-392.
- Gratton E, Fantini S, Franceschini MA, Mantulin W, Barbieri B, inventors; Photosensor with multiple light sources. US patent 5,497,769. March 12, 1996.
- 32. Stankovic MR, Maulik D, Rosenfeld W, et al. Role of frequency domain optical spectroscopy in the detection of neonatal brain hemorrhage-A newborn piglet study. J Matern Fetal Med. 2000;9:142-149.
- Suzuki S, Takasaki S, Ozaki T, Kobayashi Y. Tissue oxygenation monitor using NIR spatially resolved spectroscopy. *Proc Soc Photo Opt Instrum Eng.* 1999;3597: 582-592.
- 34. Chance B, Nioka S, Kent J et al. Time-resolved spectroscopy of hemoglobin and myoglobin in resting and ischemic muscle. *Anal Biochem*. 1988;174:698-707.

- 35. Delpy DT. Quantitation of pathlength in optical spectroscopy. *Adv Exp Med Biol.* 1989;248:41-46.
- Flock ST, Patterson MS, Wilson BC, Wyman DR. Monte Carlo modeling of light propagation in highly scattering tissue—I: Model predictions and comparison with diffusion theory. *IEEE Trans Biomed Eng.* 1989;36:1162-1168.
- Fantini S, Hueber D, Franceschini MA, et al. Non-invasive optical monitoring of the newborn piglet brain using continuous-wave and frequency-domain spectroscopy. *Phys Med Biol.* 1999;44:1543-1563.
- Thavasothy M, Broadhead M, Elwell C, Peters M, Smith M. A comparison of cerebral oxygenation as measured by the NIRO 300 and the INVOS 5100 Near-Infrared spectrophotometers. *Anaesthesia*. 2002;57:999-1006.
- Kim M, Ward D, Cartwright C, Kolano J, Chlebowski S, Henson I. Estimation of jugular venous saturation from cerebral oximetry of arterial O2 saturation during isocapnic hypoxia. J Clin Monit. 2001;16:191-199.
- 40. Schell RM, Cole DJ. Cerebral monitoring: jugular venous oximetry. *Anesth Analg.* 2000;90:559-566.
- Goldman S. Interventions based on cerebral oximetry reduce the incidence of prolonged ventilation and hospital stay in cardiac surgery patients. In: Outcomes 2004: The Key West Meeting; May 19-24, 2004; Key West, FL.
- 42. Goldman SM, Sutter F, Ferdinand F. Optimizing intraoperative cerebral oxygen delivery using noninvasive cerebral oximetry decreases the incidence of stroke for cardiac surgical patients. *Heart Surg Forum*. 2004;7:E376-E381.
- 43. Murkin JM. Perioperative detection of brain oxygenation and clinical outcomes in cardiac surgery. *Semin Cardiothorac Vasc Anesth.* 2004;8:13-14.
- Murkin JM, Adams SJ, Novick RJ, et al. Brain oxygenation in diabetic patients during coronary surgery: a randomized prospective blinded study. *Anesth Analg.* 2005; 100:SCA101.
- Moraca R, Lin E, Holmes JH IV, et al. Impaired baseline regional cerebral perfusion in patients referred for coronary artery bypass. *J Thorac Cardiovasc Surg.* 2006;131: 540-546.
- Nakamura Y, Kawachi K, Imagawa H, et al. The prevalence and severity of cerebralvascular disease in patients undergoing cardiovascular surgery. *Ann Thorac Cardiovasc Surg.* 2004:10:81-84.
- 47. Hall RA. Brain SPECT imaging and neuropsychological testing in coronary artery bypass patients: single photon emission computed tomography. *Ann Thorac Surg.* 1999; 68:2082-2088.
- Reents W, Fordyce DJ, Lee ME, et al. Cerebral oxygen saturation assessed by near-infrared spectroscopy during coronary artery bypass grafting and early postoperative cognitive function. *Ann Thorac Surg.* 2002;74:109-114.
- 49. Oki A, Ohtake H, Okada Y, Kawada T, Takaba T. Simultaneous monitoring of somatosensory evoked potentials and regional cerebral oxygen saturation combined with serial measurement of plasma levels of cerebral specific

proteins for the early diagnosis of postoperative brain damage in cardiovascular surgery. J Artif Organs. 2004;7:13-18.

- 50. Yoshitani K, Kawaguchi M, Iwata M, et al. Comparison of changes in jugular venous bulb oxygen saturation and cerebral oxygen saturation during variations of haemoglobin concentration under propofol and sevoflurane anaesthesia. *Br J Anaesth.* 2005;94:341-346.
- Schwarz G. Cerebral oximetry in dead subjects. J Neurosurg Anesthesiol. 1996;8:189-193. J Neurosurg Anesthesiol. 1997;9:76-77.
- 52. Beese U, Langer H, Lang W, Dinkel M. Comparison of near-infrared spectroscopy and somatosensory evoked potentials for the detection of cerebral ischemia during carotid endarterectomy. *Stroke*. 1998;29:2032-2037.
- Büchner K, Meixensberger J, Dings J, Roosen K. Nearinfrared spectroscopy—not useful to monitor cerebral oxygenation after severe brain injury. *Zentralbl Neurochir*. 2000;61:69-73.
- Mille T. Near infrared spectroscopy monitoring during carotid endarterectomy: which threshold value is critical? *Eur J Vasc Endovasc Surg.* 2004;27:646-650.
- 55. www.casmed.com. Accessed July 27, 2007.
- Paul B, Bo C. http://crisp.cit.nih.gov. Accessed July 27, 2007.
- 57. Strangman. Non-invasive neuroimaging using nearinfrared light. Soc of Biol Psych. 2002;52:679-693.
- McCormick PW, Stewart M, Goetting MG, Dujovny M, Lewis G, Ausman JI. Noninvasive cerebral optical spectroscopy for monitoring cerebral oxygen delivery and hemodynamics. *Crit Care Med.* 1991;19:89-97.
- 59. Ito H, Kanno I, Fukuda H. Human cerebral circulation: positron emission tomography studies [review]. *Ann Nucl Med.* 2005;19:65-74..
- 60. MacLeod DB, Ikeda K, Keifer JC, Moretti E, Ames W. Validation of the CAS adult cerebral oximeter during hypoxia in healthy volunteers. *Anesth Analg.* 2006;102:S162.
- 61. MacLeod DB, Ikeda K, Moretti E, Keifer JC, Grocott H. Using the CAS cerebral oximeter to estimate cerebral venous oxygen saturation. In: American Society of Anesthesiologists; October 22, 2005; Atlanta, GA.
- 62. Macmillan CSA, Andrews PJD. Cerebrovenous oxygen saturation monitoring: practical considerations and clinical relevance. *Intensive Care Med.* 2000;26:1028-1036.
- Chieregato A, Calzolari F, Trasforini G, Targa L, Latronico N. Normal jugular bulb oxygen saturation. J Neurol Neurosurg Psychiatry. 2003;74:784-786.
- 64. MacLeod K, Ikeda K, White W. Pilot study of fore-sight cerebral oximeter in cardiac patients. *Anesth Analg.* 2006;102:S162.
- 65. Fischer GW, Reich DL, Plestis KA, Griepp RB. The application of absolute cerebral oximetry during aortic surgery. *Anesthesiology.* 2006;105:A430.
- 66. Rodriguez RA, Cornel G, Semelhago L. Cerebral effects in superior vena caval cannula obstruction: the role of brain monitoring. *Ann Thorac Surg.* 1997;64:1820-1822.

- 67. Sung-Hee H, Chong-Sung K. Obstruction of the superior vena cava cannula detected by desaturation of the cerebral oximeter. *J Cardiothorac Vasc Anesth.* 2005;19:420-421.
- 68. Yeh T Jr, Austin EH, Sehic A. Rapid recognition and treatment of cerebral air embolism: the role of neuromonitoring. *J Thorac Cardiovasc Surg.* 2003;126:589-591.
- 69. Orihashi K, Sueda T, Okada K, Imai K. Malposition of selective cerebral perfusion catheter is not a rare event. *Eur J Cardiothorac Surg.* 2005;27:644-648.
- Hagl C, Khaladj N, Karck M, Kallenbach K. Hypothermic circulatory arrest during ascending and aortic arch surgery: the theoretical impact of different cerebral perfusion techniques and other methods of cerebral protection. *Eur J Cardiothorac Surg.* 2003;24:371-378.
- Hogue CW Jr, Palin CA, Arrowsmith JE. Cardiopulmonary bypass management and neurologic outcomes: an evidencebased appraisal of current practices. *Anesth Analg.* 2006; 103:21-37.
- 72. Yamamoto Y, Georgiadis AL, Chang H-M, Caplan LR. Posterior cerebral artery territory infarcts in the New England medical center posterior circulation registry. *Arch Neurol.* 1999;56:824-832.
- Rankin JM, Silbert PL, Yadava OP, Hankey GJ, Stewart-Wynne EG. Mechanism of stroke complicating cardiopulmonary bypass surgery. *Aust N Z J Med.* 1994;24:154-160.
- Paciaroni M, Silvestrelli G, Caso V. Neurovascular territory involved in different etiological subtypes of ischemic stroke in the Perugia stroke registry. *Eur J Neurol.* 2003; 10:361-365.
- 75. Gottesman RF, Sherman PM, Grega MA, et al. Watershed stroke after cardiac surgery—diagnosis, etiology, and outcome, *Stroke*. 2006;37:2306-2311.
- Murkin JM, Adams SJ, Novick RJ, et al. Monitoring brain oxygen saturation during coronary bypass surgery: a randomized, prospective study. *Anesth Analg.* 2007; 104:51-58.
- 77. Orihashi K, Sueda T, Okada K, Imai K. Near-infrared spectroscopy for monitoring cerebral ischemia during selective cerebral perfusion. *Eur J Cardiothorac Surg.* 2004;26:907-911.
- Casati A, Fanelli G, Pietropaoli P, et al. Continuous monitoring of cerebral oxygen saturation in elderly patients undergoing major abdominal surgery minimizes brain exposure to potential hypoxia. *Anesth Analg.* 2005;101:740-747.
- 79. Fischer GW, Reich D, Plestis KA, Griepp RB. Results utilizing absolute cerebral oximetry monitoring suggests the need for tailored patient management during cardiac surgery. In: Outcomes 2006: The Key West Meeting; May 17-20, 2006; Key West, FL.
- Reents W, Muellges W, Franke D, Babin-Ebell J, Elert O. Cerebral oxygen saturation assessed by near-infrared spectroscopy during coronary artery bypass grafting and early postoperative cognitive function. *Ann Thorac Surg.* 2002;74:109-141.